



# Turnkey solutions for surgical and telehealth

## Envoy EHR workstation with medical-grade Teguar® All-in-One



### For Surgical:

Easy to clean, no risk for interference or shock hazards



### For Telehealth:

Integrated webcam with built-in microphone & speakers



### Medical Grade Quality

Designed specifically for nurses in a surgical environment



### Dependable Battery Life

Durable, solid state drive lowers power consumption



### Fanless with IP65 Rating

Ingress protection and fanless design simplifies infection control



### 60601 Patient Safety


Complies with the highest standards of patient safety

**enovate**  
MEDICAL  
An Ergotron Company


Enovate® Envoy equipped with 22" All-in-one



## UNV00064 - 22" Widescreen



## UNV00062 - 24" Widescreen



	22" Widescreen	24" Widescreen	
<b>COMPUTER</b>	<b>Processor</b>	Intel Quad Core i5	Intel Quad Core i5
	<b>System Memory</b>	2x DDR4 SO-DIMM, 16GB	2x DDR4 SO-DIMM, 16GB
	<b>System Type</b>	64-bit Compatible	64-bit Compatible
	<b>Chipset</b>	Tiger Lake	Tiger Lake
	<b>OS Compatible</b>	Windows 11 (OS Not Included)	Windows 11 (OS Not Included)
	<b>Storage Disk Drive</b>	256GB Solid-State Drive	256GB Solid-State Drive
	<b>Security</b>	TPM 2.0	TPM 2.0
	<b>Wireless Communication</b>	Wifi 802.11 a/b/g/n/ac, Bluetooth 5.2	Wifi 802.11 a/b/g/n/ac, Bluetooth 5.2
	<b>Integrated Features</b>	2x 2W Speakers, Reading Light	2x 3W Speakers
	<b>Power Requirement</b>	DC 19 Volts 2.5 X 5.5 Connector	DC 19 Volts 2.5 X 5.5 Locking Connector
	<b>Weight</b>	11 lbs without stand	16 lbs without stand
	<b>Warranty</b>	3 Years Standard, Priority Repair*	3 Years Standard, Priority Repair*
<b>Mount</b>	75mm or 10mm VESA	75mm or 10mm VESA	
<b>DISPLAY</b>	<b>Screen Size</b>	21.5" (16:9)	23.8" (16:9)
	<b>Resolution</b>	1920 x 1080 (Full HD)	1920 x 1080 (Full HD)
	<b>Brightness</b>	250 nits	250 nits
	<b>Contrast Ratio</b>	1000:1	1000:1
	<b>Touch Screen</b>	Projected Capacitive (Multitouch)	Projected Capacitive (Multitouch)
<b>INPUT/OUTPUT</b>	<b>USB</b>	4x USB 2.0, 4x USB 3.0	4x USB 2.0, 4x USB 3.0
	<b>Ethernet</b>	2x RJ45 GbE LAN	2x RJ45 GbE LAN
	<b>Camera</b>	2 Megapixel Webcam	5 Megapixel Webcam
	<b>Audio</b>	1x 3.5mm Headphone Jack	1x 3.5mm Headphone Jack
	<b>Microphone</b>	2x Built-in	2x Built-in
	<b>Video Outputs</b>	1x HDMI	1x HDMI
	<b>Video Inputs</b>	1x HDMI	1x HDMI

\*Upon receipt at repair depot, 3-day priority repair with free return shipping

Discover the advantages of Enovate Medical's AIO Workstations for your clinical workflows.  
 Call 888-909-8920 or email us at [info@enovatemedical.com](mailto:info@enovatemedical.com)

